

ALOS



ALOS Eyelid Implants



GOLD									
Reference	TH 06	TH 08	TH 10	TH 12	TH 14	TH 16	TH 18		
Weight (gr)	0.6	0.8	1.0	1.2	1.4	1.6	1.8		
Lenght (mm)	10.5	12.5	15.5	17.5	19.5	21.5	24.5		
Width (mm)	5.7	5.7	5.7	5.7	5.7	5.7	5.7		
Thickness (mm)	0.6	0.6	0.6	0.6	0.6	0.6	0.6		
Sterilization method	Steam								
Expiration	3 years								

C€ 0318

PLATINUM										
Reference	PTH 06	PTH 08	PTH 10	PTH 12	PTH 14	PTH 16	PTH 18			
Weight (gr)	0.6	0.8	1.0	1.2	1.4	1.6	1.8			
Lenght (mm)	10.2	12.4	15.0	17.0	18.8	21.0	24.0			
Width (mm)	5.7	5.7	5.7	5.7	5.7	5.7	5.7			
Thickness (mm)	0.6	0.6	0.6	0.6	0.6	0.6	0.6			
Sterilization method	Non sterile, sterilize by moist heat at 121 °C 15 min.									

C€ 0434 REV:00

