




ALOS



The future of vision



ALOS Eyelid Implants

| Brand name | ALOS GOLD-PLATINUM EYELID IMPLANTS | | | | | | |
|----------------------|---|-------|-------|-------|-------|-------|-------|
| Indication | Eyelid implants are designed for logophthalmos functional defect treatment. Logophthalmos is defined as the inability to close the eyelids completely, both in spontaneous blinking and in forced closure. Specially indicated for patients with paralytic logophthalmos and logophthalmos that present normal passive closure. | | | | | | |
| Characteristics | <ul style="list-style-type: none">• Eyelid implants has a spherical design suitable for the curvature of the eyeball.• High biocompatible materials.• Thinner design to improve eyelid aesthetic results.• Rounded and smooth edges to avoid long term extrusion.• Suture holes with fixation channels.• 99.9% purity.• 7 standard sizes.• Different sizes or customized implants available upon request.• GOLD: Sterile / PLATINUM: Non sterile. | | | | | | |
| GOLD | | | | | | | |
| Reference | TH 06 | TH 08 | TH 10 | TH 12 | TH 14 | TH 16 | TH 18 |
| Weight (gr) | 0.6 | 0.8 | 1.0 | 1.2 | 1.4 | 1.6 | 1.8 |
| Lenght (mm) | 10.5 | 12.5 | 15.5 | 17.5 | 19.5 | 21.5 | 24.5 |
| Width (mm) | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 |
| Thickness (mm) | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 |
| Sterilization method | Steam | | | | | | |
| Expiration | 3 years | | | | | | |

CE 0318

| PLATINUM | | | | | | | |
|----------------------|--|--------|--------|--------|--------|--------|--------|
| Reference | PTH 06 | PTH 08 | PTH 10 | PTH 12 | PTH 14 | PTH 16 | PTH 18 |
| Weight (gr) | 0.6 | 0.8 | 1.0 | 1.2 | 1.4 | 1.6 | 1.8 |
| Length (mm) | 10.2 | 12.4 | 15.0 | 17.0 | 18.8 | 21.0 | 24.0 |
| Width (mm) | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 | 5.7 |
| Thickness (mm) | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 | 0.6 |
| Sterilization method | Non sterile, sterilize by moist heat at 121 °C 15 min. | | | | | | |

CE 0434

REV:00