




The future of vision



AJL VISC Sodium Hyaluronate

Reference	Viscosity (1.5 s - 125°C)	Back Stop
AJL VISC 1.4%	15000-20000 mPa.s	Blue
AJL VISC 2.0%	40000-60000 mPa.s	Pink
AJL VISC 3.0%	160000-200000 mPa.s	Green

AJL VISC is a non absorbable, highly purified and non-inflammatory viscoelastic solution of high molecular weight Sodium Hyaluronate, clear, isotonic, sterile and non-pyrogenic for intraocular injection during surgery of the anterior segment of the eye.

- High endothelium protection.
- Excellent elasticity, cohesion and coatibility properties.
- Easy to remove from anterior chamber.
- Isotonic, steril, non-pyrogenic, non-antigenic and iso osmotic viscoelastic solution
- Outstanding rheological properties.

pH	7.0-7.5
Osmolarity	250-350 mOsmol/l
Volume	1 ml
Supply	Prefilled glass syringe and disposable cannula.

CE 0434